Exploring the selectivity of C18 phases with Phenyl and PFP functionality

Introduction

HPLC columns packed with C18 bonded phases dominate all other types of bonded phases used for reversed-phase separations, and with good reason. C18 phases offer good retention and selectivity for a wide variety of sample types and have proved to be very rugged and reliable. However, the popularity of C18 bonded phases has obscured the advantages that other bonded phases, such as CN, phenyl, fluorinated, or polar embedded phases can offer. This is unfortunate because there are numerous cases where choosing a phase other than a typical C18 can lead to a more optimum separation better resolution, faster analysis, higher sensitivity, etc.

A solution to this dilemma is to add selectivity to a C18 bonded phase. This approach has been successfully applied in the development and production of C18 phases with Phenyl and pentaflurophenyl groups attached to the alkyl chain. This paper will discuss the mechanisms of separation provided by these enhanced selectivity phases and suggest how they may be used in developing reversed-phase separations.

: The effect of N, α and k on resolution (R_s)



Increasing N, α or k increases Resolution (R_s). However, as can be seen from these plots, increasing either N or k suffers from quickly diminishing returns. Increases α (selectivity), on the other hand, does not have this problem and, therefore, becomes the most powerful of these three variable to optimize when developing a separation.

: Hydrophobic Binding Interaction **FIGUR**



The dominant separation mechanism in reversed phase chromatography is hydrophobic binding interaction between solute molecules in the mobile phase and the stationary phase, i.e. the bonded phase, such as C18, C8, etc. Although the actual retention mechanism is not well understood, it has been found useful to describe the process as solute partitioning from the mobile phase into the bonded phase.







When pentafluorophenyl (PFP) groups are added to the stationary phase, an additional mechanism of separation, hydrogen bonding interaction, often plays a role in separations. In this separation of hydroxybenzoic acid isomers, hydrogen bonding interactions appear to be responsible for the peak elution order reversal on the ACE C18-PFP compared to the phases that lack hydrogen bonding interactions, i.e. the ACE C18 and ACE C18-AR.

Robert T. Moody, Thomas J. Waeghe, Carl L. Zimmerman MAC-MOD Analytical, Inc., Chadds Ford, PA

Sample Identities:	
1.6 mm 1. 4-hydroxybenzoi	c ac
H_2PO_4 , 2.3-hydroxybenzoid	c ac
3. benzoic acid	

Both Genistein and Apigenin are better retained on the ACE C18-AR and the ACE C18-PFP than the ACE C18 because of π - π interactions. However, significantly better selectivity, and thus Resolution, is achieved on the ACE C18-PFP due to shape selectivity. The greater rigidity of the C18-PFP phase creates conditions where Apigenin interacts more strongly with the stationary phase than Genisteir due the difference in their shape in solution.

FIGURE 10: Taking Advantage of Selectivity Differences in **Reversed Phase HPLC**

Popular base deactivated C18 columns are unable to provide adequate separation for all peaks in this sample of substituted methoxybenzene isomers. The additional selectivity provided by the pentafluorophenyl group in the C18-PFP phase is able to accomplish baseline separation of all 8 peaks.

Conclusion

The C18 phase remains the most popular phase for reversed phase HPLC separations because of its applicability to a wide range of sample types and because of its ruggedness and reliability. However, the lone mechanism of separation provided by a C18 phase, hydrophobic binding interactions, is sometimes not enough to achieve an acceptable separation. This poster has shown how other mechanisms of separation, i.e., π - π , hydrogen bonding, dipole-dipole, and shape selectivity, can be useful in separating peaks not well separated by hydrophobic binding interactions alone. A new generation of stationary phases that combine the benefits of C18 with the additional mechanisms of separation offered by phenyl and pentafluorophenyl can be powerful tools to use when developing new HPLC methods or improving existing methods.