

Conditions

Column: ACE UltraCore 2.5 SuperC18
 Dimensions: 150 x 3.0 mm; 2 x 150 x 3.0 mm (coupled); 3 x 150 x 3.0 mm (coupled)
 Part Number: CORE-25A-1503U
 Mobile Phase: A: 0.1% formic acid in H₂O
 B: 0.1% formic acid in MeCN

Time (mins)				
150 x 3.0 mm	300 x 3.0 mm	450 x 3.0 mm	%B	
-	0.00	0.00	5	
0.00	0.47	0.94	5	
20.00	40.47	60.94	100	
25.00	45.47	75.94	100	
26.00	46.47	76.94	5	
46.00	86.47	136.94	5	

Flow Rate: 0.43 mL/min
 Temperature: 80 °C
 Detection: UV, 254 nm

