

Conditions

Column:	ACE 3 C18
Dimensions:	50 x 2.1 mm
Part Number:	ACE-111-0502
Mobile Phase:	MeOH/0.5 mM ammonium acetate in H ₂ O (50:50 v/v)
Flow Rate:	0.2 mL/min
Injection:	20 µL
Temperature:	Ambient
Detection:	ESI MS/MS (+ve mode)
Sample:	Metabolites derivatised with 2-nitrobenzaldehyde to form nitrophenyl derivatives, prior to LC-MS analysis

